
食物有害物質的快速生物檢測方法 Fast Biosensor for Detection of Harmful Substance in Food

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
A rapid, sensitive and low cost technology for food and drug screening

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Badan-labeled biosensor in water. Left: no antibiotics added; right: penicillin G added, which shows an increase in intensity

標記了螢光素的生物傳感器。 左:未有加入抗生素:右:加入 了青霉素G,螢光明顯增強 Fluorescein-labeled biosensor in water. Left: no antibiotics added; right: penicillin G added, which shows a prominent increase in fluorescence intensity

Concern about traces of beta-lactam antibiotics (e.g. penicillins and cephalosporins) in food and in the environment is spurring the search for simple, specific, sensitive, and rapid methods of detection. We have designed a biosensor for detecting antibiotics based on enzymes (beta-lactamases) secreted by pathogenic bacteria to destroy antibiotics. The biosensor works on the principle of gross conformational changes in the enzyme when it 'sees' antibiotic substrate. We attached an environment-sensitive fluorescent molecule to a flexible loop close to the enzyme's active site. When the enzyme binds a beta-lactam antibiotic, the loop moves, changing the environment of the probe and causing the probe to fluoresce strongly. Our technology allows the fine-tuning of fluorescence color from the biosensor for different purposes. The sensor responds to beta-lactam antibiotics and beta-lactamase inhibitors. It can detect penicillin G at nanomolar levels. It can also be used for screening enzyme inhibitors and hence finds its application in drug screening.

Principal Investigator

Prof. K. Y. Wong

Department of Applied Biology & Chemical Technology
Contact Details

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
Patent Application No: 10/401,867 (USA), PCTON 2004000253 (Patent Cooperation Treaty)

Special Features and Advantages

- "turn-on" fluorescent biosensor
- sensitive: can detect antibiotic concentration as low as 10 nM in water
- rapid and easy to operate
- specific to b-lactam antibiotics and b-lactamase inhibitors
- can be adapted to "high-throughput" screening, i.e. detection of large number of samples at the same time, thus reducing time, cost and manpower
- no expensive machines are required
- low production cost

Applications

- Detection of b-lactam antibiotics in food and in the environment.
- To test whether milk or other dairy products are containing trace levels of antibiotics in order to meet the FDA requirement.
- To test the antibiotic level during the processes for making yogurt (thick fermented liquor made from milk). If the antibiotic level is too high, bacteria cannot grow and yogurt cannot be produced.
- To detect antibiotics in blood samples
- For high-throughput drug screening to discover new antibiotic components (enzyme inhibitors)

Award

Gold Award, the 5th China International Invention Expo (September 2004)